



	★ GOOD ★	★★ BETTER ★★	★★★ BEST ★★★	★★★★ BARIATRIC ★★★★★
<b>Item #</b>	P503	P703	P903	P750
<b>Product Name</b>	Prime Care Bed Model P503	Prime Care Bed Model P703	Prime Care Bed Model P903	Prime Care Bed Model P750
<b>Height Range</b>	7.8" - 30"	7" - 30"	7" - 30"	11" - 30"
<b>Weight Capacity</b>	450 lbs	550 lbs	600 lbs	750 lbs
<b>Resident Weight Capacity</b>	425 lbs	500 lbs	550 lbs	700 lbs
<b>Bed weight</b>	215 lbs	215 lbs	222 lbs	282 lbs
<b>Deck Width</b>	35"	35"	36" - 39" - 42"	39" - 42" - 48"
<b>Width option</b>	None	Opt Integrated 39", 42"	Integrated Push Button	Integrated Push Button
<b>Length</b>	76" - 80"	76" - 80" - 84"	80" - 84" - 88"	80" - 84" - 88"
<b>Length option</b>	--	Option 88"	STD Integrated	STD Integrated
<b>Mobility Heights</b>	Any	Any	Any	Any
<b>One Touch Floor Lock</b>	--	✓	✓	✓
<b>Frame Deck Design:</b>	Grid	Slat	Slat	Slat
<b>Vertical Lift</b>	✓	✓	✓	✓
<b>Motors</b>	3	4	4	4
<b>Trend /Rev-Trend</b>	--	✓ - only w/ staff	✓ - only w/ staff	✓ - only w/ staff
<b>Staff Control</b>	--	Optional	Optional	Optional
<b>Auto Contour</b>	✓	✓	✓	✓
<b>Leg/Heel Lift</b>	--	✓	✓	--
<b>Comfort (Cardiac)Chair</b>	--	✓ - only w/Staff	✓	✓
<b>Floor Lock Casters</b>	✓	✓	✓	✓
<b>All Wheels Swivel</b>	✓	✓	✓	✓
<b>Back Lite Pendent</b>	--	--	✓	✓
<b>Under Bed Light</b>	--	--	✓	✓
<b>Battery Back up</b>	--	Optional	Optional	Optional
<b>IV Pole</b>	--	✓	✓	✓
<b>Warranty</b>	<b>Frame</b>	10 years	15 years	15 years
	<b>Welds</b>	15 years	Lifetime	Lifetime
	<b>Electrical</b>	2 years	5 years	5 years

LTC Bed Accessories		
Assist Bars		
PLTCAB	PLTCAR	PLTCDR
Assist Bar - P903/P703/P503/P750	Assist Rail - P903/P703/P503/P750	Deluxe Assist Rail - P903/P703/P750

LTC Bed Accessories			
P703 Accessories			
P703WE	P703LE	P703SC04	P703-BBK
Width Extension 39"-42"	Length Extension - 88"	Staff Control	Battery Backup w/Bracket

LTC Bed Accessories				
Accessories				
P903WE	P903SC03	P903SC05	P903-BBK	P903TRP
Width Extension 39"-42" for P903/P750	Staff control - Linak motors for P903	Staff Control - Ti-Motion motors for P903/P750	Battery Backup w/ Bracket for P903/P750	Transporter for P903/P703/P503/P902/P601/P301

LTC Bed Accessories											
Bed Ends											
35BWBE-A	35BWBE-B	35BWBE-FM	42BWBE-A	42BWBE-B	42BWBE-FM	35MCBE-A	35MCBE-B	35MCBE-FM	42MCBE-A	42MCBE-B	42MCBE-FM
35" Brownlee Woods Bed End OSLO Cherry	35" Brownlee Woods Bed End OSLO Oak	P35" Brownlee Woods Bed End OSLO Figured Mahogany	42" Brownlee Woods Bed End OSLO Cherry	42" Brownlee Woods Bed End OSLO Oak	42" Brownlee Woods Bed End OSLO Figured Mahogany	35" Mill Creek Bed End OSLO Cherry	35" Mill Creek Bed End OSLO Oak	P35" Mill Creek Bed End OSLO Figured Mahogany	42" Mill Creek Bed End OSLO Cherry	42" Mill Creek Bed End OSLO Oak	42" Mill Creek Bed End OSLO Figured Mahogany