



NITRILE AND LATEX GLOVES

9+1

Amazing

PURPLE

- 92885 (XS)
- 92886 (S)
- 92887 (M)
- 92888 (L)
- 92889 (XL)



\$XX, XX
\$XX, XX NET

Sonic

COBALT BLUE

- 93775 (XS)
- 93776 (S)
- 93777 (M)
- 93778 (L)
- 93779 (XL)



\$XX, XX
\$XX, XX NET

Transform

ICE BLUE

- 9889A5 (XS)
- 9889A6 (S)
- 9889A7 (M)
- 9889A8 (L)
- 9889A9 (XL)



\$XX, XX
\$XX, XX NET

Vibrant

CREAMY BEIGE

- 98225 (XS)
- 98226 (S)
- 98227 (M)
- 98228 (L)
- 98229 (XL)



\$XX, XX
\$XX, XX NET

Promotion with ALL Free Goods program : Dealers to send monthly free goods claims to Supermax, no longer than 30 days after the end of the period. Promotions can be canceled at anytime. Supermax will not pay any distributor surcharge for free goods handling.

NEW

2+1



Biodegradable paper cup

Rolled edge for greater comfort and control | Water-repellent | BPI certificate | Water-based ink
5 oz | 150 ml

- C00006BIO
- C00007BIO
- C00008BIO
- C00009BIO
- C00010BIO
- C00011BIO



20 sleeves of 50 cups 1 000 cups per case

\$XX, XX
\$XX, XX NET

Earloop surgical masks

MADE IN CANADA

OPTIMAL STRONG



3 ply, fluid-resistant, enclosed aluminum nosepiece, high breathability, latex-free

3 PLY

2+1



Compliance Standards

ASTM F2100

BLUE

- 2110 (Level 1) **\$XX, XX \$XX, XX NET**
- 2120 (Level 2) **\$XX, XX \$XX, XX NET**
- 2130 (Level 3) **\$XX, XX \$XX, XX NET**

PINK

- 2122 (Level 2) **\$XX, XX \$XX, XX NET**
- 2132 (Level 3) **\$XX, XX \$XX, XX NET**

	ASTM1	ASTM2	ASTM3
Fluid-resistant (mmHg)	80	120	160
Bacterial filtration efficiency	≥95%	≥98%	≥98%
Particulate filtration efficiency (@0.1 µm)	≥95%	≥98%	≥98%
Delta P (mm H ₂ O/cm ²)	<5,0	<6,0	<6,0
Flame spread	Class 1		

Hand sanitizer moisturizer

2 in 1 antiseptic lotion

Dermaseptyl

8060

2 oz | 60 ml

28 bottles per case

~~\$XX, XX~~ **\$XX, XX NET**

8500

16.9 oz | 500 ml

16 bottles per case

~~\$XX, XX~~ **\$XX, XX NET**



2+1

1+1

Tray cover

Manufactured from a unique heavy-weight paper stock non-sterile

TC5125

8.5" x 12.25" | 25 cm x 31 cm

~~\$XX, XX~~ **\$XX, XX NET**



Self-sealing sterilization pouches

~~\$XX, XX~~ **\$XX, XX NET**



200 per box	50 boxes per case 10 000 pouches per case	2.5" x 4" 6 x 10 cm	SP22540
	20 boxes per case 4000 per case	2.75" x 9" 7 x 23 cm	SP27590
		5.5" x 6.5" 14 x 17 cm	SP52565
	15 boxes per case 3000 per case	3.5" x 9" 9 x 23 cm	SP35900
	10 boxes per case 2000 per case	3.5" x 5.25" 9 x 13 cm	SP35525
5.25" x 10" 13 x 25 cm		SP52510	
7" x 13" 18 x 33 cm		SP71300	
100 per box	10 boxes per case 1000 per case	10" x 15" 26 x 38 cm	SP10150
		12" x 15" 31 x 38 cm	SP12150
	8 boxes per case 800 per case	12" x 18" 31 x 46 cm	SP12180

Disinfectant with peroxide

Eliminates 99.99% of germs
Made with 0.5% (w/w) hydrogen peroxide
Fragrance-free

Per-Oxy Activated

Disinfectant for hard surfaces

SD1000

32 oz | 946 ml — DIN 02504863

12 bottles per case

~~\$XX, XX~~ **\$XX, XX NET**

Disinfecting wipe

100% biodegradable cotton

WD1000

160 wipes — DIN 02519275

8 bottles per case

~~\$XX, XX~~ **\$XX, XX NET**



PurCotton biodegradable sponges

100% cotton, 4 ply, non-woven,
lint-free, non-sterile

1122

2" x 2" | 5 cm x 5 cm

~~\$XX, XX~~ **\$XX, XX NET**

1144

4" x 4" | 10 cm x 10 cm

~~\$XX, XX~~ **\$XX, XX NET**

