

# Quarterly Promotions!

## Universal Kit

### FPSK04 shown:

- 2 handles, 1 IPR gauge
- 4 single sided (1 of each .08 mm, .10 mm, .13 mm, .18 mm)
- 6 double sided (1 of each .11 mm, .15 mm, .21 mm, .30 mm, .40 mm, .50 mm)

## Your Complete Solution for Interproximal Reduction (IPR) and Composites!

- Primary Function: IPR when .40mm and .50mm are needed, yet can still be used to contour composites
- IPR gauges for exact reduction
- Attachable handle for greater ease of use and patient comfort



FitStrip  
Instructional  
Video



## IPR Kit

### FPSK05

- 2 handles, 1 IPR gauge
- 10 double sided (2 of each .11 mm, .21 mm, .30 mm, .40 mm, .50 mm)

## For the Dentist looking to solely use for IPR.

- Double sided finishing strip ranging from .11mm to .50mm.
- IPR Gauges for exact reduction
- Attachable handle for greater ease of use and patient comfort.






**2 Free Saws**  
with purchase of a **FPSK04**  
or **FPSK05** starter kit.








2 FPSSAW-1  
value \$68

**FitStrips** are a complete line of interproximal trimmers and finishers from the company that invented tight interproximal contacts, **Garrison® Dental Solutions.**

**Buy 4  
packs  
Get 1  
FREE**

Refills	Single Sided	Double Sided
Super fine 15 Micron 	FPSXFSS-4 0.08 mm thick	FPSXFDS-4 0.11 mm thick
Fine 30 Micron 	FPSFSS-4 0.10 mm thick	FPSFDS-4 0.15 mm thick
Medium 46 Micron 	FPSMSS-4 0.13 mm thick	FPSMDS-4 0.21 mm thick
Coarse 90 Micron 	FPSCSS-4 0.18 mm thick	FPSCDS-4 0.30 mm thick
Serrated 	FPSSAW-4 0.05 mm thick	

Refills	Single Sided	Double Sided
Medium 46 Micron 	N/A	FPSTMS-4 0.40 mm thick
Medium 46 Micron 	N/A	FPSXTMS-4 0.50 mm thick
Super fine 15 Micron 	FPSGXFSS-4 0.08 mm thick	N/A
Fine 30 Micron 	FPSGFSS-4 0.10 mm thick	N/A
Medium 46 Micron 	FPSGMSS-4 0.13 mm thick	N/A

\*To redeem goods, email a scanned copy of your dealers invoice dated April 1 - June 30/2025 to: Garrison Dental Solutions, Attn: Liz Kammers to ekammers@garrisdental.com with INCA25Q2 in the Subject line. Allow 4-6 weeks for delivery.

**Garrison**

Toll-free 888.437.0032  
gds@garrisdental.com  
www.garrisdental.com

**ORDER TODAY! 888.437.0032**  
or contact your Dealer Representative.

INCA25Q2