

DiaPex® Plus

Calcium Hydroxide with Iodoform

**BUY 1 REGULAR KIT,
GET 1 REFILL!**



Pre-mixed • Highly Radiopaque • Antibacterial

[A1001-501] Regular Kit 1 x 2g syringe and 20 tips
[A1001-504] Refill 1 x 2g syringe



NEW!

DIA-X SIL™

Hydrophilic VPS Impression Material

BUY 2, GET 1 FREE!

- 1 Excellent elastic recovery
- 2 Hydrophilic properties for excellent flow
- 3 High tear strength from deformation
- 4 Long-lasting dimensional stability
- 5 Excellent detail reproduction
- 6 Patient comfort: easy to remove

[A2007-1003] Heavy Body (2 x 50 ml)

Heavy Body

[A2008-1003] Light Body (2 x 50 ml)

Light Body



DUO-GUN

Cordless Backfill Obturation Device

BUY & GET:
1 BOTTLE GP OBTURATORS
1 DIA-ROOT BIO SEALER



- Variable temperature settings: 140°C, 170°C, 200°C
- Temperature reaches 200°C within 15 seconds
- Designed for both left-handed and right-handed
- Operate hand trigger with minimal effort
- Strong and durable tips are available in both 23G & 25G



[4007-2001] Duo-Gun Regular Kit
Handpiece & Charger, GP Obturation Tips (4 x 23G), Plunger Set, Multi-Tool, Safety Caps (2), Plunger O-Ring Assembly (2), Scrubbing Brush, Adapter & Cord.

DIA ULTRASONIC

Ultrasonic Endo Activator

BUY & GET:
2 DIA-ROOT BIO SEALER KITS

- 1 Dual-frequency: 30 KHz and 45 KHz
- 2 Flexible and durable NiTi tip
- 3 Easy to use and compact
- 4 Fast charging: 3.5 hours for a full charge
- 5 Long lasting: up to 500 cycles
- 6 Powerful and efficient irrigation

[4009-1001] Regular Kit
Handpiece, charger, tips (2/box), wrench, adapter & cord



Dia-ProISO .04 & .06 Taper GP (mm-Marked)

Complements: .04/.06 Constant Tapered Files, ProFile® ISO Rotary Files, EndoSequence® Files, K3® Files, EdgeSequel™ Sapphire, EdgeSequel™ Utopia, EdgeFile® X7, Vortex®, Vortex Blue™, VFN VaryFlex® NEO, VFU VaryFlex® Universal, VFF VaryFlex® Flow



**BUY 3, GET 1 FREE!
BUY 6, GET 3 FREE!***

**of the same value and series*

.04 Taper/Size #10-80 | #ML114-600 series .06 Taper/Size #15-80 | #ML116-600 series

Dia-ISO GT .04 & .06 Taper GP



.04 Taper/Size #10-80 | #114S600 series

.06 Taper/Size #15-80 | #116S600 series

D.LUX+

Dual-Wavelength Cordless Curing Light

BUY & GET:
20 DIAFIL FREE! (6 x A1, 8 x A2, 6 x A3)

- 1 Powerful 2,400 mW/cm² output
- 2 Dual wavelength: 385 - 515 nm
- 3 360° swivel light probe
- 4 5 versatile curing modes
- 5 Built-in Light Intensity Sensor
- 6 Wireless inductive charger

[4008-1001] D-Lux+ Regular Kit

Handpiece & charger, light probe, protector, disposable sheaths (200/box), power supply (adapter)

