



## Dia-X™ Files Heat-Treated NiTi Rotary Files



\*Dia-ProT GP Assorted F1-F5 (#ML150-S693)

**BUY 6, GET 1 FREE!**  
**BUY 11, GET 2 MORE + 1 DIA-PROT GP\* FREE!**



ProTaper®  
(Equivalent Size)

- Convex triangular x-section reduces rotational friction
- Progressive taper design enhances flexibility and efficiency
- Navigate calcified or severely curved canals with control
- High corrosion resistance with coated surface
- 11 mm length handle offers improved accessibility to teeth
- Perfect match with Dia-ProT gutta percha & paper point

### Item Codes (4/pack)

- #571-100 #DX 19mm
- #571-101 #D1 21mm
- #571-102 #D2 21mm
- #571-103 #D3 21mm
- #571-104 #D4 21mm
- #571-105 #D5 21mm
- #571-191 #D1-D4 21mm
- #571-201 #D1 25mm
- #571-202 #D2 25mm
- #571-203 #D3 25mm
- #571-204 #D4 25mm
- #571-205 #D5 25mm
- #571-291 #D1-D4 25mm

## NiTi Hand Files (Nickel Titanium)



**BUY 4, GET 1 FREE!**


- Color coded handles with pre-inserted endo stops
- Produces smooth and progressive enlargement of canal
- NiTi has memory so no need to pre-bend files
- NiTi resists breakage, decreases zipping and stripping
- Increased angular rotation reduces fracture
- 6 files/box. Eligible series: **NT501, NT502, NT503**
- **K-Files:** available in 21 & 25 mm. Sizes #8-80 and ass't
- **H-Files:** available in 21 & 25 mm. Sizes #8-80 and ass't
- **Reamers:** available in 21 & 25 mm. Sizes #8-80 and ass't

## Dia-Pro W™ GOLD Gutta Percha & Paper Point

**BUY 4, GET 1 FREE!**

\*of the same value

- Precisely match **WaveOne® GOLD Rotary Files**
- Millimeter marked: quick and easy method for depth measurements
- Cadmium free, latex free, and highly radiopaque
- Length = 28 mm +/- 1
- Convenient one-hand slide open packaging



### Gutta Percha (60/box)

- ML168-S601 Small
- ML168-S602 Primary
- ML168-S603 Medium
- ML168-S604 Large
- ML168-S691 Assorted

### Paper Point (100/box)

- MP268-S601 Small
- MP268-S602 Primary
- MP268-S603 Medium
- MP268-S604 Large
- MP268-S691 Assorted

## DIA-ROOT™ Bio Sealer BioCeramic Root Canal Sealer

**BUY 1, GET 1 DIA-ISOGT GP\* FREE!**

\*Dia-ISOGT .04 #15 (114-S603)



- BioCeramic: excellent biocompatibility
- pH 12: highly antibacterial
- Superior bonding to dentin & GP
- Low solubility: zero shrinkage
- Washout resistant
- Optimal flowability of 22mm
- Saves time: pre-mixed and pre-filled


### #1003-701 Dia-Root Bio Sealer

Contains: 1 x 2g syringe, 20 tips, 1 cap

## DiaEtch™ 37% Phosphoric Acid Etching Gel

**BUY 3, GET 1 FREE!**

\*of the same value



- Fast! 15-second etch
- Precise! 25 gauge pre-bent tips
- Washes off easily with no residue
- Beautiful blue color
- Ideal consistency and biocompatible

### #A2001-3103 Regular Pack

Contains: 2 x 3 ml syr. and 10 disposable tips

### #A2001-3202 Refill Pack

Contains: 1 x 5 ml syr. and 5 disposable tips