




# Buy 4 refill packs, Get a 5th pack FREE!

## Rally™ Mini Polishers for Composite




### Remove - Choose the Blue

Polishers for fast removal and contouring



### Smooth - Choose the Red

Polishers for smoothing and silky shine



### Shine - Choose the Gray

Polishers for ultra high shine

Item	Refills	Qty.	Price
FPC010	Blue coarse grit point composite polisher	30	\$87
FPC040	Blue coarse grit cup composite polisher	30	\$87
FPCM010	Red medium grit point composite polisher	30	\$87
FPCM040	Red medium grit cup composite polisher	30	\$87
FPCF010	Gray fine grit point composite polisher	15	\$87
FPCF040	Gray fine grit cup composite polisher	15	\$87

\*Canadian MSRP

Slow speed right angle latch shank. Use with water.



Do not exceed maximum speed



Wear breathing protection



Wear eye protection



Use waterspray minimum 50ml/min



Kit	Qty.	Price
FPC01 Composite polisher complete kit	60 (10 each)	\$229

**Garrison**

order toll free 888.437.0032 | [www.garrisondental.com](http://www.garrisondental.com)  
offer expires 24/11/2024 | Cannot be combined with other advertised specials

PMCA1024OH

# Buy 4 refill packs, Get a 5th pack FREE!

Free pack must be of equal value.

## Adjust and Polish Zirconia and Lithium Disilicate the eZr™ Way!

### eZr™ Intra-oral Adjustment:



- **No chatter, chipping or microfracture** - Unlike diamonds bonded directly to a hard metal shaft, eZr's rubber substrate cushions the diamonds to prevent damage to the restoration.
- **No heat build-up** - The carefully engineered substrate does not transfer heat during use.
- **FAST!** - eZr's coarse grit blue grinders make your adjustments fast and convenient.

Item	Description	Qty.	Price
FPZC030	 Blue coarse grit disc Gross removal	5/pk	\$238
FPZC010	 Blue coarse grit flat end taper Gross removal	5/pk	\$180
FPZM010 (not in kit)	 Green medium grit small point polisher Smoothing and pre-polishing	5/pk	\$75
FPZM020	 Green medium grit large point polisher Smoothing and pre-polishing	5/pk	\$113
FPZM030	 Green medium grit disc polisher Smoothing and pre-polishing	5/pk	\$113
FPZM040	 Green medium grit cup polisher Smoothing and pre-polishing	5/pk	\$113
FPZF010 (not in kit)	 Orange fine grit small point polisher High gloss polishing	5/pk	\$75
FPZF020	 Orange fine grit large point polisher High gloss polishing	5/pk	\$113
FPZF030	 Orange fine grit disc polisher High gloss polishing	5/pk	\$113
FPZF040	 Orange fine grit cup polisher High gloss polishing	5/pk	\$113

\*Canadian MSRP



**1** Start by gross adjusting and contouring using the blue coarse grit grinders. Use the flat end taper for occlusal/labial/lingual/buccal surfaces and the disc for interproximal surfaces.



**2** Use the green medium grit polishers for smoothing and pre-polishing. Use the point for the occlusal grooves, the cup for labial/lingual/buccal, and the disc for interproximal surfaces.



**3** Use the orange fine grit polishers for high gloss polishing. Use the point for the occlusal grooves, the cup for labial/lingual/buccal, and the disc for interproximal surfaces.



Beautiful results!

Buy an eZr kit (FPZK01) Get a fine grit small point polisher (FPZF010) 5 pack

**FREE!** (\$75 value)



Kit	Qty.	Price
FPZK01 Complete system kit	2 grinders 6 polishers	\$302

Slow speed right angle latch shank Steam autoclaveable

# Garrison®

systems are trusted by  
more clinicians worldwide.



Restorative Dentistry  
Matrix Materials &  
Systems



## *Class IIs next to the canine are easier with Strata-G™*

The short ring shape enables stable placement on the mesial of the first premolar, which is usually a challenging situation for sectional matrix systems.



Scan the QR code for  
Dr. Eduard Epure full  
case study.



Case courtesy of Dr. Eduard Epure, Montreal, QC

# Buy 4 rings (2 packs of 2 rings), Get 1 ring FREE!


Free ring must be of equal value.

## STRATA-G™

### Nickel titanium strength

Nitinol drawn-wire nickel titanium produces optimal tooth separation for **tight** contacts.

### Improved anatomy – No ring spring-off



faster, more  
predictable  
flash control



Item	Qty.	Reg.
<b>SG400 Blue (universal) Ring</b> Used for Premolars, pedo; any tooth with a short crown	2/pk	\$369
<b>SG500 Orange (tall) Ring</b> Used for typical molars, taller teeth, or stacking over blue for MODs	2/pk	\$369
<b>SG600 Green (wide prep) Ring</b> Used for extremely wide preps, missing cusps	2/pk	\$369

## Composi-Tight® 3D Fusion™



**Hugs the curves** - Soft-Face™ silicone tips with enhanced marginal ridge anatomy hug the matrix band to the tooth for flash free restorations.



### Dentistry's first Wide Prep ring

gives you excellent contact and contour on the most challenging restorations.

**Get a grip!** All Composi-Tight® 3D Fusion™ separator rings incorporate Ultra-Grip™ retention extensions that anchor them securely to virtually any tooth. These rings stay where you place them!



Item	Qty.	Reg.
<b>FX400 Blue (universal) Ring</b> Used for premolars, pedo; any tooth with a short crown	2/pk	\$342
<b>FX500 Orange (tall) Ring</b> Used for typical molars, taller teeth, or stacking over blue for MODs	2/pk	\$342
<b>FX600 Green (wide prep) Ring</b> Used for extremely wide preps, missing cusps	2/pk	\$342

\*Canadian MSRP

For more styles and sizes, check out our website at [garrisdental.com](http://garrisdental.com), or call your preferred dealer.

# Buy 4 packs of Garrison® bands, Get a 5th pack FREE!

Free pack must be of equal value.

	Pediatric	Premolars	Premolars with deep cervical preps	Premolars and small molars	Molars	Molars with deep cervical preps	Item	Qty.	Price	
SG Series	Strata-G™ matrix bands with optimized marginal ridge									
	thickness 0.0016"			NA				SG050	100	\$190
	height >>	SG050 3.5 mm	SG100 4.5 mm	NA	SG175 5.5 mm	SG200 6.5 mm	SG300 7.5 mm	SG100	100	\$190
								SG175	100	\$190
								SG200	100	\$190
								SG300	100	\$190
								SG300	100	\$190

	Pediatric	Premolars	Premolars with deep cervical preps	Premolars and small molars	Molars	Molars with deep cervical preps	Item	Qty.	Price	
FXH Series	Firm matrix bands with eZ-Place™ technology									
	thickness 0.0015"	NA						FXH100	100	\$173
	height >>	NA	FXH100 4.3 mm	FXH150 4.3 mm w/ext 5.9 mm	FXH175 5.3 mm	FXH200 6.1 mm	FXH300 6.1 mm w/ext 7.7 mm	FXH150	100	\$173
								FXH175	100	\$173
								FXH200	100	\$173
								FXH300	100	\$173
								FXH300	100	\$173

	Pediatric	Premolars	Premolars with deep cervical preps	Premolars and small molars	Molars	Molars with deep cervical preps	Item	Qty.	Price	
FX Series	Full Curve Bands with non-stick coating and occlusal placement tab									
	thickness 0.0016"	NA						FX100	100	\$184
	height >>	NA	FX100 4.0 mm w/ext 4.4 mm	FX150 4.1 mm w/ext 6.0 mm	FX175 5.0 mm w/ext 5.6 mm	FX200 5.4 mm w/ext 6.6 mm	FX300 6.0 mm w/ext 8.7 mm	FX150	60	\$184
								FX175	100	\$184
								FX200	100	\$184
								FX300	60	\$184
								FX300	60	\$184

	Pediatric	Premolars	Premolars with deep cervical preps	Premolars and small molars	Molars	Molars with deep cervical preps	Item	Qty.	Price	
SM Series	Non-stick coating									
	thickness 0.0016"							SM050	60	\$175
	height >>	SM050 3.2 mm w/ext 4.2 mm	SM100 4.6 mm	SM150 3.8 mm w/ext 5.4 mm	SM175 5.5 mm	SM200 6.4 mm	SM300 6.4 mm w/ext 8.6 mm	SM100	100	\$175
								SM150	60	\$175
								SM175	100	\$175
								SM200	100	\$175
								SM300	60	\$175

	Pediatric	Premolars	Premolars with deep cervical preps	Premolars and small molars	Molars	Molars with deep cervical preps	Item	Qty.	Price	
M Series	Original bands									
	thickness 0.0015"							M050	50	\$129
	height >>	M050 3.2 mm w/ext 4.2 mm	M100 4.6 mm	M150 3.8 mm w/ext 5.4 mm	M175 5.5 mm	M200 6.4 mm	M300 6.4 mm w/ext 8.6 mm	M100	100	\$129
								M150	50	\$129
								M175	100	\$129
								M200	100	\$129
								M300	50	\$129

\*Canadian MSRP

For more styles and sizes, check out our website at [garrisdental.com](http://garrisdental.com), or call your Garrison® representative at 888.437.0032.


To redeem free goods, simply email a scanned copy of your dealers' invoice dated 1/10/24 - 24/11/24 to Garrison® Dental Solutions, Attn: Liz at [ekammers@garrisdental.com](mailto:ekammers@garrisdental.com) with **PMCA1024OH** in the subject line. Allow 4-6 weeks for delivery. **Please indicate free item code on the dealer invoice when submitting redemption.**

# Buy 4 packs of Garrison® wedges, Get a 5th pack FREE!

Free pack must be of equal value.

## Strata-G™ Wedge


A better fit by design. The hollow underside sits lower and seals deeper interproximally while the soft fins allow Strata-G® Wedges to adapt to concavities and irregularities.



Refills	Item	Qty.	Price
X-Small Yellow	SGYL	100	\$85
Small Blue	SGBL	100	\$85
Medium Orange	SGOR	100	\$85
Large Green	SGGR	100	\$85

## 3D Fusion™ Wedge


Combining Soft-Face™ adaptive materials, a firm inner core and advanced mechanical features to produce a wedge that truly works. Every time.



Refills	Item	Qty.	Price
X-Small Yellow	FXYL	100	\$83
Small Blue	FXBL	100	\$83
Medium Orange	FXOR	100	\$83
Large Green	FXGR	100	\$83

## A+Wedge® Astringent Coated

Controls gingival bleeding with an aluminum sulfate coating that has an astringent effect on the gingiva, improves clarity of restoration area and eliminates need to apply an astringent separately



Refills	Item	Qty.	Price
X-Small Yellow	GWAYL	100	\$86
Small Blue	GWABL	100	\$86
Medium Orange	GWAOR	100	\$86
Large Green	GWAGR	100	\$86

## WedgeWands®


WedgeWands® are designed to twist off the wand, but many prefer to leave it attached to help retract the cheek and for easy removal after the procedure is complete.



Refills	Item	Qty.	Price
X-Small Yellow	SWYL	100	\$74
Small Blue	SWBL	100	\$74
Medium Orange	SWOR	100	\$74
Large Green	SWGR	100	\$74

## G-Wedge™

If you prefer traditional wedge placement with either a hemostat or cotton pliers, G-Wedges give you the same great shape as WedgeWands®, just without the applicator wand.



Refills	Item	Qty.	Price
X-Small Yellow	GSWYL	100	\$64
Small Blue	GSWBL	100	\$64
Medium Orange	GSWOR	100	\$64
Large Green	GSWGR	100	\$64

\*Canadian MSRP

*All Garrison products are covered by our six-month money-back guarantee. If you are not satisfied with your product, return the unused portion and we will issue a refund.*

**Garrison**

order toll free 888.437.0032 | [www.garrisondental.com](http://www.garrisondental.com)  
offer expires 24/11/2024 | Cannot be combined with other advertised specials

PMCA1024OH

© 2024 Garrison Dental Solutions, LLC