


# NITRILE EXAM GLOVES

REF PM6-280X

ACCELERATOR FREE | CHEMO DRUG TESTED | TEXTURED FINGERTIPS | POWDER-FREE | BLUE



Blueberry™

## ABOUT THE GLOVE

Specially formulated without the use of any chemical accelerators or sulfur, PRIMED's Blueberry™ glove protects both the wearer and patient from Type I (*natural rubber latex protein*) and Type IV hypersensitivities (*Chemical Accelerators: Thiurams, Thiazoles & Carbamates*). Advanced lightweight technology provides excellent tactile sensitivity, dexterity and strength, and the fully microtextured glove surface with enhanced finger texturing provide better grip in all conditions.

## PROTECTION STANDARDS

**Exceeds the current ASTM D6319:** Standard Specification for Nitrile Examination Gloves for Medical Application.

**Tested against chemotherapy drugs to the current ASTM D6978:** Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs.  
(List of tested chemotherapy drugs available upon request.)

**Tested to ASTM F1671 for Viral Penetration.**

## FEATURES & BENEFITS

### LIGHTWEIGHT NITRILE TECHNOLOGY

Combines strength with tactile sensitivity, allowing intricate and precise movement.

### FLEXIBLE AND COMFORTABLE FIT


Moves with the wearer to facilitate natural dexterity.

### FULLY MICRO-TEXTURED WITH ENHANCED FINGERTIP TEXTURING

Provides reliable grip when handling small instruments and objects.

### LOW DERMA - ACCELERATOR FREE NITRILE

Prevents Type IV allergies (Dermatitis) in users with sensitivities.



## AVAILABLE SIZES



## QUANTITIES

GLOVES/BOX	BOXES/CASE	GLOVES/CASE
250	10	2500



Patented Technology



PRIMED Medical Products Inc.  
200, 2003-91 Street SW, Edmonton, AB, Canada, T6X 0W8  
TF +1.877.877.4633 | F +1.780.497.7670 | W PRIMED.com | E info@PRIMED.com

NA-REVA3