

ValueEtch

Etching Gel



Available in: Mini, Standard, and Mega size kits.

ValueEtch Etchant Gel is a 37% phosphoric acid etchant for enamel and dentin.

It's advanced formulation provides excellent thixotropic and handling characteristics, making it ideal for etching enamel in the selective etch technique or dentin in the total etch technique.

ValueEtch has a distinctive blue color for visible application.

- Excellent surface retention
- Easily rinsed without leaving residues
- Clearly visible due to its color
- Thixotropic



Available in:

ValueEtch Etching Gel Mini Kit

Includes: 3 ml syringe and 3 tips

Sku: [VMETCH3](#)



ValueEtch Etching Gel Standard Kit

Includes: 10 ml syringe and 5 tips

Sku: [VMETCH10](#)



ValueEtch Etching Gel Mega Kit

Includes: 50 ml syringe, 3 x 2 ml empty syringes, and 20 tips

Sku: [VMETCH50](#)

