



ValueSeal

Pit & Fissure Sealant



ValueSeal Pit & Fissure Sealant is a fluoride releasing light cured sealant that can be used in the prevention of caries on primary and permanent teeth.

The low viscosity allows for easy application and deep penetration.

ValueSeal is color changing so you can see where you place it - when light activated, the color changes to a standard shade for aesthetics.

- Fluoride release
- Radiopaque
- Very low viscosity allows easy application and deep penetration
- Color changing to facilitate placement and aesthetics
- Excellent surface wear resistance due to inclusion of particles under 0.4 microns
- Use on primary or permanent teeth



ValueSeal Pit & Fissure Sealant Kit
Includes: 2 x 2 g syringes and 5 tips

Sku: [VMSEAL2](#)



Available at Valuemed Professional Products Ltd.