


Osstell ISQ guides your implant treatments



A  COMPANY

6490 Hawthorne Drive, Windsor, ON
Canada, N8T 1J9
Phone: 1 800 265 6277
office.ca@wh.com | www.osstell.com

Osstell Beacon Brochure
Osstell AB Date of issue: 2022-03-25 25103-01 US



Go from Good to Great

Intuitive | Easy to use | Hassle-free

The Osstell Beacon tells you when an implant is ready for loading in a matter of seconds.



PREDICTABLE OUTCOMES

Our method eliminates chance. When measuring implant stability, you will get an objective value, making implant treatments more predictable - even for patients with risk factors.





REDUCED TREATMENT TIME

Immediate loading is becoming more and more common. Patients want results fast. With Osstell®, you only wait as long as needed.




PROVEN METHOD

Backed by more than 1300 scientific studies, our method is based on RFA to determine whether or not an implant is stable enough. The result is presented as an ISQ value of 1-100. The higher the ISQ, the more stable the implant.



The Osstell ISQ Scale

The Osstell Beacon uses proprietary and evidenced-based RFA technology, backed by more than 1300 scientific studies. Visit osstell.com/scientific-database to learn more.



OsstellConnect is a cloud-based solution offering a new way for clinicians to use data to enhance their dental implant treatments. It is the world's largest implant online portal with more than a quarter of a million datapoints. OsstellConnect is currently free of charge and accessible to Osstell Beacon and Osstell IDx.